- FR - Notice d'utilisation (p. 2)
- EN - Instructions for use (p. 26)
- IT - Istruzioni per l'uso (p. 48)

## Table des matières

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a) .... Begin therapy .....	36
b) .... Pause therapy.....	37
16. Stop therapy .....	37
a) .... Stop therapy manually .....	37
b) .... Stop therapy automatically.....	37
17. Cleaning, disinfection, storage and return .....	38
a) .... Preparations for cleaning.....	38
b) .... Cleaning and disinfection of WoundEL® device .....	38
c).... Cleaning WoundEL® master cable / dressing cable.....	39
d) .... Cleaning WoundEL® AC/DC adapter .....	39
e) .... Return of products .....	39
18. Charging the batteries.....	39
a) .... General information.....	39
b) .... Battery indication.....	39
c).... How to charge the batteries .....	40
19. View THERAPY HISTORY .....	40
a) .... View THERAPY HISTORY .....	40
b) .... Return to Home from THERAPY HISTORY .....	41
20. Set Time, Date, Patient ID and Language .....	41
a) .... Access ADVANCED SETTINGS .....	41
b) .... Set Time .....	41
c).... Set Date .....	41
d) .... Set Patient ID.....	41
e) .... Set language .....	42
f).... Return to HOME from ADVANCED SETTINGS.....	42
21. Access SERVICE for advanced settings.....	42
a) .... Access SERVICE .....	42
b) .... View SYSTEM INFORMATION .....	42
c).... View ERROR COUNTERS .....	42
d) .... Delete errors.....	43
22. Technical description .....	43
a) .... Maintenance and repair.....	43
b) .... Electrical isolation .....	43
c).... Essential performance.....	43
d) .... Error messages/technical trouble shooting .....	43
e) .... Electromagnetic compatibility .....	44
23. Symbols .....	46

## Introduction

These instructions for use are designed to help you use the WoundEL® electrostimulation device correctly. Please read the following information carefully. Keep these instructions in case you need to use them again. Should they be lost, all notices can be downloaded from the website [www.woundel.com](http://www.woundel.com) Failure to comply with the instructions for use may have serious consequences, including personal injury. In the event of a serious incident occurring in connection with the device, the manufacturer must be notified of the incident, and must notify the competent authority of the Member State in which the patient and/or user is established. If in doubt about how to use it, contact a healthcare professional. The characteristics (frequency, intensity) of the therapy will be given by the responsible clinician.

### 1. Technical specifications

#### a) WoundEL® device (ESDWO1)

<b>Environmental conditions</b>	
<b>Ambient pressure</b>	
Operation	70 kPa to 106 kPa
Transport and storage	70 kPa to 106 kPa
<b>Ambient temperature</b>	
Operation	+5 to +40°C
Transport and storage	-25 to +70°C
<b>Ambient humidity</b>	
Operation	15% to 93% non-condensing
Transport and storage	Up to 93% non-condensing
<b>Dimensions</b>	190x130x45 mm
<b>Weight</b>	0.6 kg
<b>Classification</b>	
Electrical protection	Class II
Applied parts	Type BF. Applied parts include dressing and disperser electrodes, master cable plugs and sockets, dressing cable plugs and sockets, clips. Power supply unit not to be subjected to the requirements of applied parts.
Protection against harmful ingress	IP22
Integrated software protection	Software not connected, USB port
Mode of operation	Continuous
Electromagnetic compatibility (Emission)	Class B
Service life	5 years
CE mark according to MDD (93/42/CEE)	CE 0459
Device classification according to MDD	Ila

#### b) AC/DC adapter / battery charger

Model / Type	DTF n°YEWBC1
Input	100-240 Vac, 50 - 60 Hz, 0.6 A
Output	18 V DC, 1.0 A, 18 W
IP class	IP 52
Double insulation Class II	
<b>Environmental conditions</b>	
<b>Ambient temperature</b>	
Operation	0° to +40°C
Transport and storage	-25°C to +70°C
<b>Humidity</b>	0% to 90% non-condensing
The expected service life of the battery charger is 2 years.	

#### c) Internal battery

Type	Lithium-ion
Voltage	14.4 V (4 cells of 3.6 V each)
Capacity	2600 mAh

**d) Patient current pulse**

The values below are only valid with use together with the WoundEL® dressing and disperser electrodes.

Pulse type	Monophasic direct current (DC)
Pulse frequency	Adjustable, 128 Hz or 64 Hz ± 20 %
Pulse duration	Non adjustable, 140 microseconds ± 20 %
Polarity	Adjustable, negative or positive
Intensity	Adjustable, 0 to 42 milliampere (mA) ± 20 %
Therapy duration	Non adjustable, preset to max 30 min
Load impedances range	80 to 1000 Ω

**2. System description**

These Instructions for Use are intended to give technical instructions to the users in handling the device. The dosage and therapy instructions shall be given by the responsible clinician.

**a) System overview**

To connect the system, see Paragraph 12. "Preparations for use".


Figure 4. System overview

**b) Pack content**


Figure 5. WoundEL® device


Figure 6. AC/DC adapter

Figure 7. Master cable


Figure 8. Dressing cable with clip for WoundEL® dressing electrode (2 cables supplied)


Figure 9. Storage/transport bag

**c) Disposables (not supplied with the device)**


Figure 10. WoundEL® dressing electrode


Figure 11. WoundEL® disperser electrode

#### d) Spare parts

Accessories	Specifications	Part number
Dressing electrode 15x15 cm		EWOPAN1
Disperser electrode		EWODIS1
Master cable	Length: 180 cm (+3/-6)	YEWOMC1
Dressing cable with clip		YEWODCE2
AC/DC adapter/battery charger		YEWOBC1
Storage/transport bag		YEWOTB1

### 3. Destination, target population and users

#### a) Purpose of the device

The WoundEL<sup>®</sup> device is a medical device for the treatment of chronic wounds using electrostimulation.

#### b) Target population

The target population is any person over the age of 18 who requires electrostimulation treatment to help heal wounds and who does not have any of the contraindications listed in chapter 7.

#### c) Users

Users are professional carers, patients or family members in acute care centres (hospitals), follow-up care centres or geriatric centres, but also at home for patients cared for in their own homes (home care).

### 4. Use

The WoundEL<sup>®</sup> system is intended to be used:

- For activating and supporting wound healing by low frequency monophasic direct current impulses generated by the electrical device.

The recommendations for using the WoundEL<sup>®</sup> system are as follows:

- Frequency of sessions: two sessions of 30 minutes per day
- Maximum duration of treatment: 30 days according to available clinical data. Use beyond 30 days must be evaluated by the practitioner through a benefit / risk ratio.
- It is up to the practitioner to define and adjust the duration of use and the various parameters of the sessions according to the progress of the treatment and its evaluation

### 5. Clinical benefit and expected performance

The clinical benefits of the WoundEL<sup>®</sup> electrostimulation system are:

- Significant pain reduction
- Significant wound size reduction

The claimed performances of the WoundEL are the following:

- Promotes granulation,
- Promotes epithelialization,
- Increases transcutaneous oxygen tension in the tissue,
- Increases capillary density in the tissue,

### 6. Indications

- Pressure ulcers, stage 3 to 4
- Leg ulcers
- Diabetic foot ulcers

### 7. Contraindications and side effects

#### a) Contraindications

- Wound entirely covered with dry necrotic tissue. Adequate debridement shall be done according to clinical practice prior to start of WoundEL<sup>®</sup> electrical stimulation treatment.
- Malignant wounds or malignancies anywhere between the dressing electrode and the disperser electrode.

- Untreated osteomyelitis in immediate wound proximity or anywhere between the dressing electrode and the disperser electrode.
- Enteric fistulas.
- Non-enteric fistulas which are not possible to debride/cleanse at dressing changes.
- Metal implants in the wound or anywhere between the dressing electrode and the disperser electrode when applied on the body.
- Patients with a cardiac pacemaker, connected to a defibrillator, a high-frequency surgical device or TENS device.
- Pregnancy.
- Untreated deep leg vein thrombosis and arterial occlusion in the immediate vicinity of the wound or anywhere between the dressing electrode and the disperser electrode.
- Patients on treatment with thrombolytics.
- Wounds above the mammillary line and directly over the heart.
- Do not use on patients with known allergy/sensitivity to the WoundEL® dressing electrode or its components.
- Children under 18 years.
- Patients with haemophilia.
- Patients with a known history of Parkinson's disease, epilepsy or organ fistulas.

#### b) Possible side effects

- Sensitivity/allergy to the WoundEL® dressing electrode or its components.
- Maceration of the wound edges/surrounding skin.
- Redness of the skin beneath the WoundEL® electrodes.
- Neurological sensations such as „tingling“ in the area of the WoundEL® dressing electrode.
- In very rare cases, there may be an increased sensation of pain.
- Excessive granulation:  
During the electrical stimulation treatment period, monitor the wound closely at every dressing change for signs of excessive granulation (hypergranulation). If this occurs, stop the electrical stimulation treatment and consult a health care professional for adequate treatment. Once the excessive granulation is treated, the electrical stimulation may be started again depending on the wound condition.
- Bleeding: WoundEL® is an active and effective method of wound treatment which, among other things, increases blood circulation in the wound area and surrounding skin. The increased occurrence of bleeding during electrical stimulation or during a change of dressing is usually a sign that the treatment is starting to have an effect. If the bleeding is unusually frequent or intense, please consult the treating doctor immediately. Also remember to consult *Paragraph 6 "Contraindications"* of these Instructions for Use.
- The prescriber must inform the patient of the residual risks, contraindications and possible side effects described in this leaflet.

## 8. Warnings, precautions and safety instructions

#### a) Warnings

For assistance, if needed, in setting up, using or maintaining the WoundEL device and components; or to report unexpected operations or events, please ask your local distributor or directly WoundEL® Health Care.

**Any serious incident occurring in connection with the device must be notified to the manufacturer and to the competent authority of the Member State in which the user/patient is established.**

#### General

- The WoundEL® device must only be used according to the prescriptions by the responsible clinician.
- The WoundEL® device must only be used in accordance with the information provided in these Instructions for Use.
- Do not touch the conductive part of the two electrodes and the patient simultaneously.
- Do not touch the cable connectors and the patient simultaneously.
- No modification of this equipment is allowed.
- Repair and/or modification of the WoundEL® device or any accessory by anyone other than qualified service personnel may significantly compromise the unit's ability to conduct therapy and/or void the equipment warranty.
- The settings must always be adjusted according to the instruction by the responsible clinician.
- If the device is intended to be powered by batteries, the device must be charged prior the therapy.
- Replacement of the batteries by inadequately trained personnel could result in a hazard (such as excessive temperatures, fire or explosion).

- The WoundEL® device must not be connected to a PC or any other equipment than those specified in paragraph 3. "System description".
- If the device is to be used by multiple patients the device must be disinfected according to the instructions, see *Paragraph 16 «Cleaning, disinfection, storage and return»*.
- Do not use the electrodes at the head (including the mouth), directly on the eyes, on the front of the neck, on the thorax or upper back. The electrodes must not be positioned so that the current passes through the heart.

#### **Usage environment**

- Unless the device is used within specified environmental conditions and locations, there may be a risk that the device does not operate according to specifications.
- To avoid the risk of fire or explosion, do not use the device in the presence of flammable anesthetics and/or oxygen rich environment.
- The use of accessories and cables other than those specified, with the exception of cables sold by the manufacturer of the WoundEL® device, as replacement parts for internal components, may cause incorrect function and risk of affecting other equipment due to increased emissions or decreased immunity.
- Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in these Instructions for Use, see *Paragraph 21 "Technical description"*.
- Portable and mobile RF communications equipment can affect medical electrical equipment.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations and walkie talkies can affect this equipment and should be kept at least a distance away from the equipment in accordance with *Paragraph 21.d « Electromagnetic compatibility »*.
- Placing the device in direct sunlight may cause a risk over overheating.

#### **Interactions with other treatment devices**

- Simultaneous connection of a patient to a high frequency surgical equipment may result in burns at the site of the WoundEL® electrodes and possible damage to the WoundEL® device.
- Operation in close proximity (e.g. 1 m) to a shortwave or microwave therapy medical electrical equipment may produce instability in the WoundEL® device output.
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- Particular care should be exercised generally when using any radiotherapy equipment at the same time. WoundEL® device should be observed to verify normal operation.

#### **Interactions with imaging devices**

- Because of the particular composition of the WoundEL® dressing electrode and the WoundEL® disperser electrode, false X-ray, computerized tomography (CT) or magnetic resonance imaging (MRI) images may be obtained. The WoundEL® dressing and disperser electrodes should be removed before the use of radiological diagnostic equipment i.e x-ray and replaced for the duration of use by an alternative, X-ray permeable wound covering.

#### **b) Precautions**

- There are at present no scientifically documented findings on the use of WoundEL® on children and infants, on patients with haemophilia, on patients with a known history of Parkinson's disease, epilepsy or in organ fistulas.
- Do not use on patients with known sensitivity to the dressing electrode or its components.
- To avoid the risk of infection while using the WoundEL® system, the user must maintain good cleaning practice according to cleaning and disinfection instructions set out in this manual.
- Any electrodes that have current densities exceeding 2 mA/cm<sup>2</sup> may require the special attention of the operator.
- When the device and its accessories are taken out of use, make sure that local regulations on product disposal are considered and followed. This includes electronics, batteries, cables, plastics and contaminated waste.
- Patients treated with antiplatelet agents or anticoagulants have a low risk of complications, but should be monitored in accordance with current clinical guidelines.

#### **c) Safety instructions**

- Ensure that the user, operator and responsible organization have received sufficient training on the device and on wound management by electrostimulation specifically. User training shall be performed regularly as needed and in line with the user's basic knowledge.

- Prior to allowing patients to use the device without supervision, patient assessment must be conducted by a health care professional. The assessment shall include an evaluation confirming that the patient is capable of reading and understanding the information in this instruction and on the display.
- The WoundEL<sup>®</sup> device must only be used during the environmental conditions and in the locations specified in this manual.
- Ensure that the device is properly assembled prior to using the device.
- Ensure that the cables are positioned away from the patient and operator so that the risk of strangulation and further suffocation is minimized.
- Ensure to handle the cable clips and other small parts so that the risk of swallowing, inhalation and suffocation is minimized.
- Ensure that markings on the device and in the display are readable, otherwise do not use the device.
- Do not position the WoundEL<sup>®</sup> device so that it is difficult to operate the disconnection device, i.e the power supply unit and its mains plug.
- Use only the original WoundEL<sup>®</sup> components (see *Paragraph 9 « WoundEL<sup>®</sup> dressing electrode » and Paragraph 10 « WoundEL<sup>®</sup> disperser electrode »*). These are designed to work together in such a way as to provide the maximum possible safety and efficacy.
- Never immerse the WoundEL<sup>®</sup> components in liquid.
- Do not operate the WoundEL<sup>®</sup> device if it is damaged.
- If fault occurs during operation, switch the WoundEL<sup>®</sup> device off, disconnect the cables and inform the responsible care provider.
- Follow the cleaning and disinfection instructions carefully, see *Paragraph 16 « Cleaning, disinfection, storage and return »*. Failure to follow these instructions might cause infections.
- Ensure the correct positioning and contact of the WoundEL<sup>®</sup> components and connections during therapy.
- Ensure that the correct electrode cable is connected to corresponding outlet on the master cable connection port. Incorrectly connected cables may result in the therapy not being delivered as intended.
- Children may change settings and disconnect the cables causing therapy to stop or be incorrect. Keep the WoundEL<sup>®</sup> components out of the reach of children during both therapy and storage of the device.
- The treatment frequency should be reduced to 64 Hz for wounds that are considerably smaller than the active surface area of the WoundEL<sup>®</sup> dressing electrode.
- Check the electrodes and cables for degradation. Degradation can cause therapy to not be delivered as intended.

## 9. WoundEL® device description

The WoundEL® device delivers an adjustable electrical current distributed through the cables and the wound dressing to the wound. The electrical current returns to the device via the disperser electrode and thus closes the electrical loop.


Figure 12. The WoundEL® device


Figure 13. The display

## 10. WoundEL® dressing electrode

**CAUTION** The WoundEL® must always and without exception be used with the WoundEL® dressing electrode.

WoundEL® dressing electrode is intended to administrate WoundEL® electro stimulation currents to the body from the electro stimulation device. Additionally the dressing electrode provides a physical protection of the wound.

Information about product description, indications, usage, storage and disposal can be found in the user instructions for the dressing electrode which comes with the dressing electrode products.

To connect the dressing electrode, see Paragraph 12 « Preparations for use ».


Figure 14. WoundEL® dressing electrode

## 11. WoundEL® disperser electrode

 **CAUTION** *The WoundEL® must always and without exception be used with the WoundEL® disperser electrode.*

WoundEL® disperser electrode returns WoundEL® electro stimulation currents from the body to the electro stimulation device.

Information about product description, usage, storage and disposal can be found in the user instructions for the disperser electrode which comes with the disperser electrode products.

To connect the WoundEL® disperser electrode, see Paragraph 12 « Preparations for use ».


Figure 15. WoundEL® disperser electrode

## 12. Safety check before use

- Check the WoundEL® device visually for any damage before use. If covers are damaged or any parts are loose, do not use the device.
- Check the cable and connections for pinching or damage before and after use. Pets and pests may damage the insulation of the cables. Do not use damaged cables.
- Place the device on a flat surface.
- Check safety-related labelling on the WoundEL® device for legibility.
- Ensure that markings on the device and in the display are readable, otherwise do not use the device.

## 13. Preparations for use

### a) General instructions

- Handle the WoundEL® electrostimulation system with gloves
- Always arrange the treatment cable between the WoundEL® device and the patient in such a way that no one can trip over it and so that it cannot cause any skin pressure injury to the patient. Have the device placed on a flat surface, e.g. a table, next to the patient.
- Ensure that the cables are positioned away from the patient and operator so that the risk of strangulation is minimized.
- Ensure to handle the cable clips and other small parts so that the risk of swallowing or inhalation is minimized.
- Do not position the WoundEL® device so that it is difficult to operate the disconnection device, i.e the power supply unit and its mains plug.
- Use only the original WoundEL® components (see Paragraph 9 « WoundEL® dressing electrode » and Paragraph 10 « WoundEL® disperser electrode »). These are designed to work together in such a way as to provide the maximum possible safety and efficacy.
- Only connect or disconnect any cable to or from the WoundEL® electrodes when the WoundEL device is switched off.

## b) Cables connection

1. Apply the WoundEL® dressing and disperser electrode according to the Instructions for Use, see figure 16.


Figure 16. Apply WoundEL® wound dressing and disperser electrode at a distance of 30 cm

2. Connect the master cable to the WoundEL® device: insert the round plug of the master cable into the socket, without twisting, until there is an audible click. The arrow on the cable connector must be pointing upwards, see figure 17. (To disconnect the master cable, pull the connector out without twisting.)


Figure 17. Connect the master cable to the WoundEL® device

3. Insert the green plug of the dressing cable into one of the sockets on the master cable marked with a green dot, see Figure 18. (Up to two dressing electrodes can be connected.)


Figure 18. Connect the dressing cable to the master cable

4. Connect the dressing cable to the dressing electrode using the clip, see Figure 19. The clip can be connected either direction.


Figure 19. Connect the dressing cable to the dressing electrode

5. Connect the disperser electrode to the master cable: insert the black plug of the disperser cable into the master cable socket marked with a black dot, see Figure 20.


Figure 20. Connect the disperser electrode to the master cable

 **Safety related tip:** If the cables are accidentally disconnected, there will be a message in the display and an audible alert, see 21 .c. Error messages /technical trouble shooting.

## 14. Set polarity and frequency

- The treatment frequency should be reduced to 64 Hz for wounds that are considerably smaller than the active surface area of the WoundEL® dressing electrode.
- Set the polarity and frequency according to the prescription.

### a) Access SETTINGS

1. Press ON/OFF to turn the device on.
2. At HOME (see figure 21), press  to access MENU (see figure 22).
3. Select Settings with .
4. Press  to access SETTINGS (see figure 23).


Figure 21. Home


Figure 22. Menu


Figure 23. Settings

## b) Set polarity

### Safety related tip

The polarity is preset to negative (-).

1. Select Polarity with .
2. Press . The field for polarity highlights in green.
3. The polarity can be set to either positive (+) or negative (-). Adjust the polarity with and .
4. Press . The polarity is now set.

## c) Set frequency

### Safety related tip

The frequency is preset to 128 Hz.

1. Select Frequency with .
2. Press . The field for frequency highlights in green.
3. The frequency can be set to either 64 Hz or 128 Hz. Adjust the frequency with and .
4. Press . The frequency is now set.

## d) Return to HOME from SETTINGS

1. Select Back with .
2. Press . You are back at MENU.
3. Press . You are HOME.

## 15. Start therapy

- Ensure that the user has received sufficient training on the device.
- Do not operate the WoundEL® device if it is damaged.
- If an error occurs during operation, switch the WoundEL® device off, disconnect the cables and inform the responsible health care provider.

### a) Begin therapy

1. Press ON/OFF to turn the device on, see Figure 24.

### Safety related tip

The WoundEL® device can be powered by either the mains power supply or the internal battery.

To connect the AC/DC converter, see Paragraph 17 « Charging the batteries ».


Figure 24. ON/OFF button

2. Press . The dialogue Start therapy? is shown and the field OK is highlighted, see *Figure 25*.

Figure 25. Dialogue screen Start therapy?

3. Press . The therapy is started.

4. Increase current intensity with until it starts to feel uncomfortable/painful, then decrease one (1) mA with , see *figure 26*. The intensity can be set from 1-42 mA.

Figure 26. Adjust intensity

5. The device will notify if the setting is below 5 mA.

6. The flashing circle (see *Figure 27*) in the display and the yellow LED light on the front panel (see *Figure 28*) indicate that current is being sent to the patient.


Figure 27. Flashing circle in display


Figure 28. LED indicator

### b) Pause therapy

1. To pause therapy, press . The message "paused" is shown, see *Figure 29*.

2. Resume therapy by pressing .


Figure 29. Paused

## 16. Stop therapy

### a) Stop therapy manually

1. To stop therapy manually: press . The message "Paused" is shown, see *Figure 30*.

2. Select "Stop" with .

3. Press . THERAPY INFORMATION is shown see *figure 31*.

4. To turn the device off, press . The message « The device is switching off » is shown, see *Figure 32*.

Figure 30. Paused


Figure 31. Therapy information

Figure 32. The device is switched off

### b) Stop therapy automatically

If not terminated earlier the therapy will stop automatically after 30 minutes.

THERAPY INFORMATION is shown, see *Figure 31*.

To turn the device off: press . The message « The device is switching off » is shown - see *Figure 32*

## 17. Cleaning, disinfection, storage and return

- Plastic bags, Instructions for Use, instructions, storage/transport bag and transport box must be stored outside the patient room.
- Only use cleaning and disinfecting solutions which are suitable for use on plastics and metals.
- DO NOT use chemicals or disinfectants containing alcohol, acid or alkali to clean or disinfect the WoundEL® device or the other components of the system. These can attack the materials chemically and cause technical problems and even complete failure of the device or its components.
- Always adhere to the concentrations and application times specified by the manufacturer of the cleaning and disinfecting solutions used! Please read the information provided with the packaging and the Instructions for Use. Failure to comply with these requirements can result in inadequate cleaning and disinfection of the WoundEL® components and thus to health problems for the user or patients.
- Always dry the WoundEL® device and WoundEL® cables in air at room temperature. DO NOT dry the WoundEL® components in warming cabinets after cleaning and disinfection! DO NOT autoclave! Excessive heat can cause technical problems and even complete failure of the WoundEL® components.
- Before any cleaning and disinfection, switch the electrical stimulation device off using the ON/OFF key. The display should now be completely blank. If the device is being charged, disconnect the AC/DC adapter from the mains socket and remove the AC/DC adapter cable from the electrical stimulation device. If a Master Cable is still connected to the device, always remove this from the device before cleaning and disinfection.
- DO NOT spray directly on or immerse the WoundEL® device, the WoundEL® cables or WoundEL® AC / DC adapter in water or other liquids. This can cause technical problems and even complete failure of the WoundEL® components.
- Statutory warranty rights are invalidated if the WoundEL® components have been cleaned or disinfected incorrectly or negligently and this causes technical problems. WoundEL® cannot provide any guarantee that the required level of cleaning or disinfection is achieved.
- Please always comply with the hygiene standards applicable in your medical facility with respect to cleaning and disinfection of medical devices. If in any doubt, please contact the person responsible for hygiene at your medical facility in the first instance.

### a) Preparations for cleaning

- Hygienic hand disinfection with an alcohol skin disinfectant
- Clean and disinfect the work surface.

### Protective clothing

- Protective overalls (waterproof),
- Facemask,
- Protective glasses,
- Gloves.

### b) Cleaning and disinfection of WoundEL® device

#### Preparation

1. Check the WoundEL® device on all sides and the connection sockets for visible signs of dirt see *Figure 33*.
2. Carefully detach visible dirt from the WoundEL® device and/or the connection sockets mechanically using a suitable tool (e.g. a soft brush), see *Figure 34*.


Figure 33. Check the device on all sides.


Figure 34. Detach visible dirt with a suitable tool

#### Cleaning

1. Thoroughly wipe the device with wipes soaked in disinfectant or ready to use disinfecting wipes that follow the requirements mention in this chapter, see *Figure 35*. Pay attention particularly to edges, narrow corners and ports.
2. Once the wipe is not wet anymore, throw it away and use a new one.
3. Wipe off all visible contamination/residues.


Figure 35. Wipe the device with cloth

4. Clean carefully the ports. Use a cotton-swab for the hard to access areas.
5. Once the device is visually clean, wipe it with a clean and dry swab in order to remove the remaining agent fluid and dirt.

### **Disinfection**

1. Thoroughly wipe device once again with wipes soaked in disinfectant or ready to use disinfecting wipes. Pay attention particularly to edges, narrow corners and ports. Replace the wipe when this one is not wet enough.
2. Let the device air dry during the contact time recommended by the supplier.
3. After drying the WoundEL<sup>®</sup> device, store it in the WoundEL<sup>®</sup> storage/transport bag separately packed from other WoundEL<sup>®</sup> components until next use, see Figure 36.


Figure 36. Put the device into storage/transport bag

#### **c) Cleaning WoundEL<sup>®</sup> master cable / dressing cable**

1. Check the WoundEL<sup>®</sup> cables on all sides and the connection sockets for visible signs of dirt.
2. Carefully detach visible dirt from the WoundEL<sup>®</sup> cables and/or the connection sockets mechanically using a suitable tool (e.g. a soft brush).
3. Follow the cleaning and disinfecting instructions in section 17.b.
4. After drying loosely roll up the WoundEL<sup>®</sup> cables and store it in the WoundEL<sup>®</sup> storage/transport bag separately packed from other WoundEL<sup>®</sup> components until next use.

#### **d) Cleaning WoundEL<sup>®</sup> AC/DC adapter**

1. Check the AC/DC WoundEL<sup>®</sup> adapter on all sides and the connection sockets for visible signs of dirt.
2. Carefully detach visible dirt from the AC/ DC adapter and/or the connection sockets mechanically using a suitable tool (e.g. a soft brush).
3. Follow the cleaning and disinfecting instructions in section 17.b
4. After drying the AC/ DC adapter, store it in the WoundEL<sup>®</sup> storage/transport bag separately packed from other WoundEL<sup>®</sup> components until next use.

#### **e) Return of products**

1. Ensure that the WoundEL<sup>®</sup> device and its accessories have been cleaned and disinfected according to instructions above.  
**NOTE!** Do not send any dressing or disperser electrodes when returning the WoundEL<sup>®</sup> device and its accessories.
2. Place the WoundEL<sup>®</sup> device and its accessories in the storage/transport bag and place the storage/transport bag in the cardboard box.

## **18. Charging the batteries**

- Do not charge the WoundEL<sup>®</sup> device while placed in the storage/transport bag or other types of closed or non-ventilated spaces.
- Only use the charger, i.e the AC/DC adapter delivered with the WoundEL<sup>®</sup> device.
- The AC/DC adapter and its cable should be handled carefully as to prevent damage.
- When using the AC/DC adapter, ensure that its cable is positioned between the apparatus and the socket in such a way that there is no risk of injury from tripping over it.
- If the device batteries have been fully discharged, the device menu will restart from the default settings: “Date = 01.03.2000, Daytime = xx:yy, Language = English” once the batteries are recharged.

#### **a) General information**

The WoundEL<sup>®</sup> device can be powered by either the mains power supply or the internal battery.

The WoundEL<sup>®</sup> device includes rechargeable lithium-ion batteries.

Since all batteries degenerate over time the recommendations below will ensure that the battery capacity is maximized during its lifetime.

#### **b) Battery indication**

The battery symbol in the upper right corner of the display indicates the status of the battery. If the battery needs charging the battery symbol will indicate low battery, see Figure 37. If the battery becomes critically low the display will give an error code (see Figure 38) and the device will shut down.


Figure 37. Low battery symbol.

Figure 38. Low battery error

### Safety related tip

The capacity of the batteries is not negatively affected if it is charged when partially empty and do not need to be fully charged at each occasion. Therefore it is only an advantage to charge the batteries as soon as an opportunity arises.

**If the device is not used for longer periods, the battery needs to be charged at least every 6 months.**

### c) How to charge the batteries

1. Connect the AC/DC adapter cable to the DC input connector on the device, see Figure 39. Make sure that the plug is properly inserted.
2. Connect the AC/DC adapter to mains supply plug in the wall.
3. Check that the battery charging indicator is shown in the upper right corner of the display, see Figure 40. The charging indicator will be shown whenever charging is ongoing, also when the device is switched off. A crossed over battery symbol in the upper right corner or full screen when switched off indicates that the charging process has stopped, see Figure 41. This could happen when the temperature is too high in the device in order to protect the batteries.


Figure 39. DC input


Figure 40. Battery charging indicator


Figure 41. Battery charging error

### Safety related tip

During the charging process the device and the AC/DC adapter heat up slightly. This is normal.

If the internal temperature in the battery is too high the charging process will stop (crossed over battery symbol). As soon as the temperature has decreased to an acceptable level the charging will continue.

New fully charged batteries will give an estimated time of treatment of 5 hours.

### Charging time

Charge	Charging time
To cover 2 therapy sessions	40 minutes
Fully charged batteries	3 h

## 19. View THERAPY HISTORY

### a) View THERAPY HISTORY

1. At HOME (see figure 42), press .
2. MENU (see Figure 43) is shown. Here you can adjust Settings, History and Advanced settings.
3. At MENU, select History with .
4. Press .
5. THERAPY HISTORY (see figure 44) is shown. Here you can view therapy sessions.
6. Select 1 of 20 with .
7. Press .
8. Select the number with and . (see Figure 44).


Figure 42. Home


Figure 43. Menu


Figure 44. Therapy history / therapy session

### b) Return to Home from THERAPY HISTORY

1. Press
2. Select Back with
3. Press . You are now back at the MENU.
4. Press . You are now HOME.

## 20. Set Time, Date, Patient ID and Language

### a) Access ADVANCED SETTINGS

1. At ADVANCED select Settings with and .
2. Press to access ADVANCED SETTINGS (See Figure 45). Here you can adjust the Time, Date, Patient ID and Language.


Figure 45. Advanced Settings

### b) Set Time

1. In ADVANCED SETTINGS (See Figure 45) Select Time.
2. Press . The field for hours highlights in green.
3. Adjust the hours with and .
4. Press . The field for minutes highlights in green.
5. Adjust the minutes with and .
6. Press . The time is now set.

### c) Set Date

1. In ADVANCED SETTINGS (See Figure 45) Select Date.
2. Press . The field for year highlights in green.
3. Adjust the years with and .
4. Press . The field for month highlights in green.
5. Adjust the month with and .
6. Press . The field for day highlights in green.
7. Adjust the day with and .
8. Press . The date is now set.

### d) Set Patient ID

1. In ADVANCED SETTINGS (See Figure 45) Select Patient ID.
2. Press . The first field highlights in green.
3. Adjust the letter with and .
4. Press . The second field highlights in green.
5. Adjust the number with and .
6. Press . The third field highlights in green.

- Adjust the number with  and .
- Press . The patient ID is now set.

### e) Set language

- In ADVANCED SETTINGS (See Figure 45) Select Language.
- Press . The field highlights.
- Adjust the language with  and .
- Press . The language is now set.

### f) Return to HOME from ADVANCED SETTINGS

- Select Back with .
- Press  to return to ADVANCED.
- Select Back with .
- Press  to return to MENU.
- Select Menu with .
- Press  to return to HOME.

## 21. Access SERVICE for advanced settings

### a) Access SERVICE

- At ADVANCED (see figure 46) select Service with .
- Press .


Figure 46. Advanced

- Service is shown (See Figure 47). Here you can access System information, Error count, Clear errors, Export error log and Clear error log.


Figure 47. Service

### b) View SYSTEM INFORMATION

- In SERVICE (See figure 47), select System info with .
- Press . SYSTEM INFORMATION is shown (See figure 48)
- To return to SERVICE: press 


Figure 48. System Information

### c) View ERROR COUNTERS

- In SERVICE (See figure 47), select Error Count with .
- Press . ERROR COUNT is shown (See figure 49)
- To return to SERVICE: press 


Figure 49. Error count

#### d) Delete errors

1. In SERVICE (See figure 47), select Clear Errors with . The dialog Clear errors? Is shown (See figure 50)
2. Select OK with .
3. Press . The message Errors cleared is shown (See Figure 51). The errors are now cleared.
4. To return to SERVICE: press sur .


Figure 50. Clear errors?

Figure 51. Errors cleared

## 22. Technical description

#### a) Maintenance and repair

- The WoundEL® device does not require preventive technical maintenance.
- Any repair needed shall be done by WoundEL® Health Care technical service or service providers approved by WoundEL® Health Care. Service providers shall fulfil service trainings set up by WoundEL® Health Care.

#### b) Electrical isolation

Isolation from supply mains is provided by the double insulated AC/DC adaptor/battery charger.

#### c) Essential performance

The transmitted signal must not exceed 50 mA.

#### d) Error messages/technical trouble shooting

Error message	Explanation	Cause	Proposed action
E0	Internal failure	Device failure	Restart the device and if the error persists, contact your distributor.
Failure E1	Internal failure	Device failure	Restart the device and if the error persists, contact your distributor.
Failure E2	Software update failure	Device failure or USB memory failure	Contact your distributor for a software update.
Failure E3	Internal failure or software update failure	Device failure or incomplete software update	Contact your distributor for a software update.
Error E4: System power failure. Please restart your device.	System power failure	Device failure	Restart the device and if the error persists, contact your distributor.
Error E5: Low battery. Please charge your device and restart it.	Battery low causing the device to power off	Low battery level	Connect charger to continue therapy
Error E6: Therapy hardware error. Please restart your device.	Internal failure	Device failure	Restart the device and if the error persists, contact your distributor.
Error E7: Error in stored data. Please restart your device.	Internal failure	Device failure	Restart the device and if the error persists, contact your distributor.
Error E8: General error. Please restart your device.	General errors	Device failure	Restart the device and if the error persists, contact your distributor.
Error: Self-test failed. Please restart your device.	Internal failure	Device failure	Restart the device and if the error persists, contact your distributor.

Crossed over battery symbol	Charging of batteries has stopped	Internal temperature too high or battery failure	Let the device cool down. Do not charge the batteries when placed in the storage/transport bag. If the error persists, contact your distributor.
Connection failure. Check cables and electrodes	Electrodes or cables insufficiently connected.	Partial loss of electrical connection to patient	Check cables and electrodes connection.
No connection. Connect cables and electrodes	Electrodes or cables not connected	Electrodes or cables not connected	Connect cables and electrodes
Possible short circuit. Check your cables and dressing.	Impedance in electrical circuit low	Electrode overlap or incorrect connection setup	Check your cables and dressing
The battery level is low. Connect charger to continue therapy.	Not possible to run therapy on battery	Low battery level	Connect charger to continue therapy

In case of any other problem, contact your distributor.

### e) Electromagnetic compatibility

Electro medical devices require special precautions against EMC and must be installed and set to operate in accordance with the EMC information provided in this manual.

**WARNING:** Avoid using this unit near or stacked onto other equipment as this may cause a malfunction. If this is necessary, this appliance and the other appliances should be monitored for normal operation.

**WARNING:** Use of accessories, transducers and cables other than those specified or supplied by the manufacturer of this device may cause an increase in electromagnetic emissions or a decrease in the immunity of this device and cause improper operation.

**WARNING:** Portable RF communication devices (including peripherals such as antenna cables and external antennas) should not be used closer than 12 inches (30 cm) to any part of the WoundEL device, including cables specified by the manufacturer. Otherwise, performance may be impaired.

**WARNING:** If the device malfunctions due to electromagnetic disturbances, the performance of the device may be impaired and result in delayed patient care.




Wireless communication devices such as wireless home networking devices, mobile phones, cordless phones and their bases, walkie-talkies can affect the device, and it is recommended to keep them at a distance of at least 0.5 m from the device.

Guidance and manufacturer's declaration – electromagnetic emissions	
Emissions test	Compliance
RF emissions - CISPR 11	Group 1
RF emissions - CISPR 11	Class B
Harmonic emissions - IEC 61000-3-2	Class A
Voltage fluctuations / flicker emissions - IEC 61000-3-3	Complies

Guidance and manufacturer's declaration – electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact	± 8 kV
	± 15 kV air	± 15 kV
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV
	± 1 kV for input/output lines	± 1 kV
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s)
	± 2 kV line(s) to earth	± 2 kV line(s) to earth
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U <sub>T</sub> ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0% U <sub>T</sub> ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
	0% U <sub>T</sub> ; 1 cycle AND 70% U <sub>T</sub> ; 25/30 cycles Single phase : at 0°	0% U <sub>T</sub> ; 1 cycle AND 70% U <sub>T</sub> ; 25/30 cycles Single phase : at 0°

Voltage interruptions IEC 61000-4-11	0% U <sub>T</sub> : 250/300 cycles	0% U <sub>T</sub> : 250/300 cycles
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz	10 V/m 80 MHz to 2,7 GHz
	385 MHz at 27 V/m 450 MHz at 28 V/m 710, 745, 780 MHz at 9 V/m 810, 870, 930 MHz at 28 V/m 1720, 1845, 1970 MHz at 28 V/m 2450 MHz at 28 V/m 5240, 5500, 5785 MHz at 9 V/m	385 MHz at 27 V/m 450 MHz at 28 V/m 710, 745, 780 MHz at 9 V/m 810, 870, 930 MHz at 28 V/m 1720, 1845, 1970 MHz at 28 V/m 2450 MHz at 28 V/m 5240, 5500, 5785 MHz at 9 V/m
Conducted RF IEC 61000-4-6	3V 0,15 MHz – 80 MHz 6V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% MA at 1kHz	3V 0,15 MHz – 80 MHz 6V in ISM bands and amateur radio bands between 0.15 MHz and 80 MHz 80% MA at 1kHz
Local magnetic fields IEC 61000-4-39	30 kHz at 8 A/m 134.2 kHz at 65 A/m 13.56 MHz at 7.5 A/m	30 kHz at 8 A/m 134.2 kHz at 65 A/m 13.56 MHz at 7.5 A/m

## 23. Symbols

Symbols	Definitions
	Compliance with EU Regulation 2017/745 of 5 April 2017 on medical devices, with identification of the notified body
	Manufacturer
	Medical Device
	Catalogue number
	Unique serial number « YY MMDD CCC » where « YYY » refers to year manufacturing, « MM » refers to month manufacturing, « DD » refers to day manufacturing, « CCC » refers to a consecutive number
	Unique Device Identification
	Transport and storage temperature limits
	Atmospheric pressure limitations
	Humidity limitation
	Keep dry
	Fragile : handle with care
	Follow Instructions for Use
	Type BF applied part
	For indoor use only
	Class II equipment (on the AC/DC adapter)
	USB connection (only for manufacturer's use)
<b>IP22</b>	This means that the device is protected against insertion of fingers and will not be damaged or become unsafe during a specified test in which it is exposed to vertically or nearly vertically dripping water
	When the LED below is lit, current is delivered to the patient
	WEEE symbol
 <b>WARNING</b>	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury
 <b>PRECAUTION</b>	Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury
 <b>Safety related tip</b>	Indicating useful information about the safe use of the device

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